

CHILDREN TOOTHBRUSHES



Private Label & Licensing




Examples...

All company and product names are trade names or trademarks and are the property of the respective owners.

YES - YOU CAN!

- Around half the toothbrushes we make are in our customer's own designs.
- It's Quick & Easy to do...!
- Available LED colours: Red, Yellow, Green & Blue
- MOQ is 5,000 pcs
- Lead time: 4-6 weeks ...from artwork approval



Factory Audits



Top Models

Flashing Timer

page
3



Alpha Model

Flashing Timer

page
5



Lambda Model

Flashing Timer

page
7



Gamma Model

Musical Timer

page
9



Kappa Model

Electric Toothbrush

page
11



Omega Model

Sonic Toothbrush

page
13



Beta Model

Electric Toothbrush

page
15



Epsilon Model

Contact Details

Twinklers Ltd

Building #37, Industrial Estate #3

Li-Song-Lang, Gong Ming

Shenzhen 518106, China

+86 755 2712-6212

+86 755 2712-6047

European Sales Office

Unit 3,

165 Granville Road

London

NW2 2AZ

United Kingdom

+44 20 8457-7000

www.bbrite.net

sales@twinklers.com



FLASHING TOOTHBRUSH



Alpha Model

Flashing Timer

2 Minutes (1 & 3 Minutes also available)
Indicating the correct brushing time.

Activation

Press button to start red LED flashing - auto off

Bristles

Tynex® - DuPont™
Soft - Grade: 612 Nylon level filaments
Number of Tufts: 22
End Rounded: 90%

Material

Neck: ABS920
Grip: ABS757
Button: TPE

Backing Card

400 gsm single sided (3.5g)

Blister Cover

PET - Thickness: 0.3mm (3.3g)

Batteries

2 x LR41 (AG3) Alkaline Button Cell - Mercury Free
Included & not replaceable
Duration: 4 months + ...lasts lifetime of bristles
Shelf life: 3 years

Size


Toothbrush: Actual - 145 x 24 (widest) mm (19g)
On Backing Card - 221 x 42 mm (26g)

Colours

032 C 801 C 802 C 807 C 210 C 803 C



Suction Cups Available





Actual Size

FLASHING TOOTHBRUSH

Lambda Model



Flashing Timer

2 Minutes (1 & 3 Minutes also available)
Indicating the correct brushing time.

Liquid Filled

Clear purified water - anti freeze, antimicrobial
and antifungal added.
Not harmful if ingested.

Bristles

Tynex® - DuPont™
Soft - Grade: 612 Nylon level filaments
Number of Tufts: 22

Material

Handle: ABS920
Character: PVC
End Cap (Suction Cup): TPR50
Pellets: AS

Backing Card

400 gsm single sided (3.5g)

Blister Cover

PET - Thickness: 0.3mm (3.3g)

Batteries

2 x LR41 (AG3) Alkaline Button Cell - Mercury Free
Included & not replaceable
Duration: 4 months + ...lasts lifetime of bristles
Shelf life: 3 years


Size

Toothbrush: Actual - 160 x 21 (widest) mm (31g)
On Backing Card - 221 x 42 mm (38g)

Colours

032 C 801 C 802 C 807 C 210 C 803 C

Suction Cups Available





Actual Size

FLASHING TOOTHBRUSH

Gamma Model



Flashing Timer

2 Minutes (1 & 3 Minutes also available)
Indicating the correct brushing time.

Activation

Press button to start LED flashing - auto off

Bristles

Tynex® - DuPont™
Soft - Grade: 612 Nylon level filaments
Number of Tufts: 22

Material

Handle: ABS920
Shrink Wrap: PET
End Cap (Suction Cup): TPR50

Backing Card

400 gsm single sided (3.5g)

Blister Cover

PET - Thickness: 0.3mm (3.3g)

Batteries

2 x LR41 (AG3) Alkaline Button Cell - Mercury Free
Included & not replaceable
Duration: 4 months + ...lasts lifetime of bristles
Shelf life: 3 years

Size

Toothbrush: Actual - 160 x 24 (widest) mm (20g)
On Backing Card - 221 x 42 mm (27g)

Colours

032 C 801 C 802 C 807 C 210 C 803 C

Suction Cups Available



ARTWORK

Printing Area



HANDLE



PREVIEW



BACKING
CARD

21

FLASHING
DEMO'S

22

SAFETY
COMPLIANCE

22



Actual Size

MUSICAL TOOTHBRUSH

Kappa Model



Musical Timer

2 Minutes (1 & 3 Minutes also available)
Indicating the correct brushing time.

Playback

Looped 1 minute medley - London Bridge, Old MacDonald, Farra Jacka
Custom soundtrack option available (.wav file required).

Activation

Press button to start music playing - auto off

Bristles

Tynex® - DuPont™
Soft - Grade: 612 Nylon level filaments
Number of Tufts: 22

Material

Handle: ABS757
Button: TPR50
Back Cover: ABS757

Backing Card

400 gsm single sided (3.5g)

Blister Cover

PET - Thickness: 0.3mm (3.3g)

Batteries


2 x LR44 (AG13) Alkaline Button Cell - Mercury Free
Included & not replaceable
Duration: 4 months + ...lasts lifetime of bristles
Shelf life: 3 years

Size

Toothbrush: Actual - 145 x 26 (widest) mm (24g)
On Backing Card - 221 x 42 mm (32g)

Colours

032 C 801 C 802 C 807 C 210 C 803 C





ELECTRIC TOOTHBRUSH

Omega Model



Activation

Press top button to start. Press bottom button to stop.

Bristles

Tynex® - DuPont™
Soft - Grade: 612 Nylon level filaments
Number of Tufts: 18
Non replaceable head

Material

Round Head: POM
Fixed Clamp: POM
Handle: ABS757
Button: TPE
End Cap: ABS757

Backing Card

250 gsm double sided (5.6g)

Blister Cover

PET - Thickness: 0.4mm (8.1g)

Motor

Voltage: 1.5v
Speed: 8,000 RPM

Batteries


1 x AA Battery - Mercury Free
Included & replaceable
Duration: 4 months + ...lasts lifetime of bristles
Shelf life: 3 years

Size

Toothbrush: Actual - 196 x 31 (widest) mm (83g)
On Backing Card - 270 x 56 mm (96g)

Colours

032 C 801 C 802 C 807 C 210 C 803 C





Actual Size

Beta Model



Activation

Press button to start vibration. Press button to stop.

Bristles

Tynex® - DuPont™
Soft - Grade: 612 Nylon level filaments
Number of Tufts: 22
End Rounded: 90%

Material

Handle: ABS757
Head: ABS757
Cradle: ABS757
Cover: ABS757
Base: ABS757
Button: TPE

Display Box

400 gsm single sided (11g)

Motor

Voltage: 1.5v
Speed: 11,000 RPM

Batteries


1 x AAA Alkaline Button Cell - Mercury Free
Included & replaceable
Duration: 2 months +
Shelf life: 3 years

Size

Toothbrush: Actual - 159 x 22 (widest) mm (40g)
In Display Box - 160 x 56 x 23 mm (56g)

Colours

032 C 801 C 802 C 807 C 210 C 803 C White





Actual Size

Epsilon Model



Activation

Press button to start. Press button to stop.

Bristles

Tynex® - DuPont™

Soft - Grade: 612 Nylon level filaments

Number of Tufts: 27

Non replaceable head

Material

Bristle Base POM

Head: ABS757

Clip Ring: ABS757

Handle: ABS757

Button: TPE

Battery Cover: ABS757

Display Box

400 gsm single sided (25g)

Motor

Voltage: 1.5v

Speed: 8,000 RPM

Batteries

1 x AA Battery - Mercury Free

Included & replaceable

Duration: 2 months

Shelf life: 3 years


Size

Toothbrush: Actual - 212 x 31 (widest) mm (80g)

In Display Box - 240 x 72 x 32 mm (124g)

Colours

032 C 801 C 802 C 807 C 210 C 803 C White





Actual Size

FLASHING TOOTHBRUSH

Theta Model



Flashing Timer

2 Minutes (1 & 3 Minutes also available)
Indicating the correct brushing time.

Activation

Press End Cap to start LED flashing - auto off

Bristles

Tynex® - DuPont™
Soft - Grade: 612 Nylon level filaments
Number of Tufts: 22

Material

Handle: ABS920
End Cap: TPR50

Backing Card

400 gsm single sided (3.5g)

Blister Cover

PET - Thickness: 0.3mm (3.3g)

Batteries

2 x LR41 (AG3) Alkaline Button Cell - Mercury Free
Included & not replaceable
Duration: 4 months + ...lasts lifetime of bristles
Shelf life: 3 years

Size

Toothbrush: Actual - 145 x 16 (widest) mm (14g)
On Backing Card - 221 x 42 mm (21g)



Colours

032 C 801 C 802 C 807 C 210 C 803 C

Suction Cups Available



BLISTER STICKER





Actual Size

REGULAR TOOTHBRUSH

Tau Model



Bristles

Tynex® - DuPont™

Soft - Grade: 612 Nylon level filaments

Number of Tufts: 22

Material

Handle: ABS757

Back Cover: ABS757

Backing Card

400 gsm single sided (3.5g)

Blister Cover

PET - Thickness: 0.3mm (3.3g)

Size


Toothbrush: Actual - 145 x 24 (widest) mm (11g)

On Backing Card - 221 x 42 mm (19g)

Colours



ARTWORK



HANDLE



PREVIEW



BACKING
CARD

21

SAFETY
COMPLIANCE

22

Backing Card



221mm

42mm



GB Designed to teach children how long to brush their teeth for. This timer toothbrush flashes for 2 minutes. Push the button to activate. Auto - off. Batteries included and non-replaceable. Makes brushing a fun part of the daily routine! Ergonomically designed handle for better control. End-rounded Tynex® bristles from DuPont™. Dentists recommend replacing your toothbrush every 3 months. Supervise brushing to ensure correct brushing technique. Do not chew bristles. For brushing teeth only. This toothbrush is not a toy and should always be used with adult supervision. Warning: Not suitable for children under 3 years. Contains small parts which could represent a choking hazard. Do not submerge in alcohol (e.g. Hand Sanitizer). Please retain this information for future reference.

DE Mit Freude lernen wie man richtig putzt. Diese Kinderszahnbürste blinkt für 2 Minuten, der empfohlenen Putzdauer. Knopf betätigen um zu aktivieren. Schaltet automatisch ab. Batterien sind enthalten, nicht auswechselbar. Spaß beim Zähneputzen, morgens und abends! Ergonomisch geformter Griff für optimalen Komfort. Abgerundete Tynex® Borsten von DuPont™. Zahnärzte empfehlen die Zahnbürste nach 3 Monaten Gebrauch auszutauschen. Üben Sie das tägliche Putzen zusammen mit Ihrem Kind. Nur zum Zähneputzen, nicht kauen! Nicht geeignet für Kinder unter 3 Jahren! Die Zahnbürste enthält verschluckbare Kleinteile - Erstickungsgefahr! Die Zahnbürste nicht mit Alkohol in Berührung bringen (z.B. Desinfektionsmittel). Bewahren Sie diese Information für den weiteren Gebrauch auf.

FR Conçu pour permettre aux enfants de se laver les dents le temps recommandé. Le minuteur clignote pendant 2 minutes. Poussez le bouton pour l'activer. S'arrête automatiquement. Les piles sont incluses et non remplaçables. Le Brossage des dents devient une activité ludique. Le bout des poils Premium DuPont™ Tynex® sont de forme arrondie. Le manche ergonomique convient aux mains des enfants. Pour le brossage des dents uniquement. Ne pas mordre. Les dentistes recommandent de changer de brosse à dents tous les 3 mois. Surveillez le brossage des dents afin de vous assurer de la bonne maîtrise du geste. Attention: déconseillé aux enfants de moins de 3 ans. Contient de petits éléments entraînant un risque d'étouffement. Garder cette notice d'information pour référence future. Ne pas immerger dans l'alcool.

LOT #
lift flap

b.brite
brush right™

TWINKLERS LTD.
LONDON NW2 2LD, UK
www.bbrite.net

DuPont™ and Tynex® are trademarks or registered trademarks of E.I. du Pont de Nemours and Company and its affiliates.



AG3 (2 pcs)

Directive 2006/66/EC
Waste Batteries
Remove batteries
before discarding &
dispose in battery
recycling bin



ITEM No. F-86C



α **Made in China** v10
EC Registration Design No.002647982-0001

Actual Size

Flashing Demo



12pc SHELF TRAY
WITH
FLASHING DEMO

Safety Compliance

All b brite toothbrushes are fully tested. They comply with the following safety regulations and labelling standards (where relevant):

EN ISO 20126 – Dentistry Performance

EN71 Part 1 – Physical Properties

EN71 Part 2 – Flammability

EN71 Part 3 – Toxic & Element Migration

EN62115 – Electrical Toy Safety

EMC – Electromagnetic Compatibility

IEC 60825-1 – LED Safety

PAH – Polycyclic Aromatic Hydrocarbons

EN ISO 3679 – Liquid Flammability

RoHS – Restriction of Hazardous Substances

RoHS – Restriction of Hazardous Substances

EC REG 10/2011 & 1935/2004 – Overall Migration

REACH REG (EC) 1907/2006 – Cadmium, Phthalate & Organotin

Microbiological Safety

Antimicrobial Preservation

Saliva and Perspiration Resistance

Class IPX4 Water Resistant

Colourfastness

GB 19342-2003 – Toothbrush Health & Safety (PRC)

